



TECHNICAL SPECIFICATIONS

Lifting Capacity: 200 lb (90 kg)

Vehicle Types: CUVs, hatchbacks, minivans, Sedans, SUVs, full-size vans

Rely on Bruno's Lifter to raise and stow a folding manual wheelchair, travel scooter or travel powerchair into your vehicle's trunk or hatchback. Connect the docking device, press a button, and your mobility device is lifted for you. Manually guide it into the back of your vehicle, and you're ready to head to your destination.

- Power raise/lower, manual rotation
- Custom bases specific to your vehicle
- Docking device included
- Fold-down arm when lift is not in use
- Eligible for auto manufacturer potential rebates*
- 3-year limited warranty

** Contact Bruno dealer for details.*

Learn more at www.bruno.com



MEDmobility

5711- Calgary Trail NW
Edmonton AB T6H 2K1
P: 780-437-3300
F: 780-430-9381
www.medmobility.ca



LIFTER FEATURES



Lifter folds down when not in use.

MEDmobility

5711- Calgary Trail NW
Edmonton AB T6H 2K1
P: 780-437-3300
F: 780-430-9381
www.medmobility.ca



Bruno Independent Living Aids, Inc.
1780 Executive Drive, Oconomowoc, WI 53066
www.bruno.com

© Bruno Independent Living Aids, Inc. 2015

MKT-L-12